



Testing of disinfectants according to
DIN 68 861, Part 1

QM H I
Revision: 1.06

Tested disinfectants

| Name | Concentration | Applicable |
|------------------------------|---|------------|
| Chloramin T | 2.5% aqueous solution | ✓ |
| Terralin protect | 0.5% aqueous solution and 100% solution | ✓ |
| Descosept Special | 100% solution | ✓ |
| Descosept | 100% solution | ✓ |
| Incidin Plus | 0.5% aqueous solution and 100% solution | ✓ |
| Mikrobac forte | 2 % aqueous solution and 100 % solution | ✓ * |
| Biguamed perfekt N | 5% aqueous solution and 100% solution | ✓ |
| perform (powder concentrate) | 1% and 3% aqueous solution | ✓ |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

* NOT SUITABLE for the net of the **KayserBett TOBI** and the net of the **travel cot NILS!**